



PROMISE
Preparing for RSV Immunisation
and Surveillance in Europe

PROMISE RSV-LabNet Protocol Library

PCR - RSV A/B, Influenza B/Vic and B/Yam

Updated 18 January 2023

This protocol has been made available by:

Dr. Adam Meijer and Dr. Lance Presser
National Institute of Public Health and the Environment,
Diagnostic Laboratory for Infectious Diseases and Screening
Department Virology
PO Box 1
3720 BA Bilthoven
The Netherlands
e-mail: adam.meijer@rivm.nl and lance.presser@rivm.nl

The copyright of this protocol resides with the original author(s).
The original author(s) are responsible for the correctness of the contents of this protocol.
This protocol may be used by the members of the PROMISE project in the context of the RSV Laboratory Network in Europe.
Distribution of this protocol outside the network and publication in any form are prohibited without obtaining explicit permission of the original author(s), unless the protocol has been published in a journal.

This project has received funding from the Innovative Medicines Initiative 2 Joint Undertaking under grant agreement 101034339. This Joint Undertaking receives support from the European Union's Horizon 2020 research and innovation programme and EFPIA. This document reflects the author's view and that neither IMI nor the European Union, EFPIA are responsible for any use that may be made of the information contained therein.

1. Reagents

RNA Extraction

MagNa Pure 96 System (Roche Material Number: 065 410 890 01) using MagNA Pure 96 DNA and Viral Nucleic Acid Small Volume Kit (Roche Catalog Number: 065 435 880 01)

Lightcycler PCR

PCR reactions should be performed on Lightcycler 480 with TaqMan™ Fast Virus 1-Step Master Mix (ThermoFisher Scientific Catalog Number: 444 443 2)

Primer Sets

RSV-A

2050F TGA ACA ACC CAA AAG CAT CA
2117R CCT AGG CCA GCA GCA TTG
2086probe2 5'-6Fam AAT TTC CTC ACT TCT CTA GTG TAG TAT TGG G 3'-BHQ1

RSV-B

17fw2 GAT GGC TCT TAG CAA AGT CAA GTT GA
120 TGT CAA TAT TAT CTC CTG TAC TAC GTT GAA
PB45probe 5'-Texas Red TGA TAC ATT AAA TAA GGA TCA GCT GCT GTC ATC CA
3'-BHQ1

Influenza B-HA (Specific B-Yam and B-Vic probes)

HA-444F2 ACC CTA CRA AAT TGG AAC CTC AGG
HA-444F3 ACC CTA CAG ACT TGG AAC CTC AGG
HA-524R ACR GCC CAA GCC ATT GTT G
Yam501probe 5'-Atto425 AAA TCC GAT TTT ACT GGT AG 3'-MGBBeclipse
Vic499probe 5'-Atto532 ATC CTT TTC CAT TGG TAA 3'-MGBFQ

Primers and Probes Information

RSV-infB mix	Label	Filter (Nm)	pmol in PCR-mix	Subset
RSV A 2050F			15	RSVA
RSV A 2117R			15	
RSV A 2086P2	FAM-BHQ1	483-533	5	
RSV B 17 fw2			15	RSVB
RSV B 120			15	
RSV B PB45	TXR-BHQ2	558-610	5	
INFB-HA-444F2			9	VIC
INFB-HA-444F3			9	YAM
INFB-HA-524R			12	
INFB-Yam501Pr	ATTO425-MGBBeclipse	450-500	3	YAM
INFB-Vic499Pr	ATTO532-MGBBeclipse	523-568	2	VIC

Target	Sequence	Filter (nm)	Subset name
RSVA	FAM-BHQ1	483-533	RSVA
RSVB	TXR-BHQ2	558-610	RSVB
InfB-Vic	ATTO532-MGBBeclipse	523-568	VIC
InfB-Yam	ATTO425-MGBBeclipse	450-500	YAM

2. Procedure

Preparation

1. Make appropriate dilutions of positive controls.
2. Isolate RNA by MagNA Pure 96 DNA and Viral Nucleic Acid Small Volume Kit (Roche Catalog Number: 065 435 880 01) Preheat two thermoblocks on 50°C and 95°C
3. Make PCR-mix for RSVA, RSVB, InfB-Vic, InfB-Yam (reagents lab)

PCR-Mix	conc.	[end-conc.]	µl
4x Taqman Fast Virus MM			5.0
RSV-InfB primer/probe mix			3.0
PCR Grade Water			7.0
Total volume			15.0

Sample Addition, Reverse Transcription, and PCR

1. Aliquot 15 ul portions of PCR-mix in 96-well plate according to plate layout template
2. Add 5 uL RNA to 15 ul PCR-mix
3. Seal the 96-well plate
4. Centrifuge briefly
5. Incubate directly 15 minutes at 50°C (heating block)*
6. Incubate directly 2 minutes at 95°C (heating block)*
7. Centrifuge briefly
8. Keep plate at 4°C if you can't run directly
9. Run PCR Program on LightCycler 480
10. Analyse results

***Note** – RIVM uses thermoblock reverse transcription step instead of performing this step in the LightCycler to allow waiting times between preparation of the reaction 96-well plate and putting the plate in the LightCycler. If the plate is only kept cool during waiting time this reduces the performance of the RT-PCR. If immediately after preparing the reaction plate the reverse transcription is performed and after that the plate is kept cool, there is no impact on performance. If a lab has enough capacity to put the plate after reaction preparation immediately on the LightCycler, this separate thermoblock reverse transcription step can be avoided. However, in that situation, the reverse transcription step has to be programmed in the temperature programme of the LightCycler.

PCR-program LightCycler 480:

Detection Format:	Multi Color Hydrolysis Probe				
PCR Program	Segment Number	Temp (°C)	Hold Time (sec.)	Slope (°C/sec.)	Acquisition Mode
Reverse Transcription	1	50	900	EXTERNAL	
Inactivate RT, activate Pol.	1	95	120	EXTERNAL	
Denaturation	1	95	60	4.4	None
Amplification (Cycles: 50)	1	95	10	4.4	None
	2	55	30	2.2	Single
Cooling	1	40	30	4.4	None

Positive Control Information*

Control	Name
Positive Control	Current strains from own repository

***Note** – The positive control material used at RIVM are strains from the RIVM repository. RIVM can provide Influenza B viruses, but at this moment, RIVM cannot provide RSV.

Additional Note(s)

Assay validation report is available if requested. The PCR performs with 100% correct results at RIVM in QCMD influenza, RSV, and respiratory panels, WHO EQAP influenza panels, and in WHO RSV EQA pilot panel by NEQAS.

Research Use Only