

Principles and Benefits of Citrate Anticoagulation

Buddy-Joe Paris

Clinical Nurse Specialist


Fresenius Medical care



**Edinburgh Critical Care, Research Group:
*What's new in ICU***


The Use of Citrate

- Citrate has been used for many years across a variety of medical applications, for example, it is widely used for the storage of red blood cells ²²
- As early as 1990, clinicians began to see the possibilities for utilising citrate as a regional anticoagulant in CRRT ¹
- Fresenius has extensive experience with regional citrate anticoagulation, with the commercial launch of Ci-Ca[®] in February 2006 after many years of research and development




The Citrate-Calcium Complex

Anticoagulation with citrate utilises a process called '**chelation**'




Citrate forms a complex with the Ca^{2+} ions, making them unavailable as a co-factor within the clotting cascade

Calcium and the Clotting Cascade




What Are We Doing in Citrate Anticoagulation?



Fresenius Ci-Ca[®] CVVHD Circuit

- Citrate solution is introduced after the access connection
- Calcium solution is infused just prior to the return connection
- Anticoagulation is achieved in the circuit, but the patient's clotting status is unaffected



Fluids for Ci-Ca[®] Therapy



Ci-Ca Dialysates
calcium free, 4 types

4% Sodium Citrate
Contains 136 mmol/L citrate



Calcium Chloride
100 mmol/L

Fluids for Citrate Therapy

Name	Language/Art. No. D/GB/F/I/E/P/NL/ DK/S/N/FIN/CZ/ SRB/SK/TR	Language/Art. No. RUS/EST/LT/LV/ GR/H	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ⁺⁺ mmol/L	Mg ⁺⁺ mmol/L	Cl ⁻ mmol/L	HCO ₃ ⁻ mmol/L	Glucose g/L
Ci-Ca® Dialysate K2	9689201	F00003194	133	2	0	0.75	116.5	20	1
Ci-Ca® Dialysate K4	F00000431	F00003195	133	4	0	0.75	118.5	20	1

Name	Language/Art. No. D/GB/F/I/E/P/NL/ DK/S/N/FIN/CZ/ SRB/SK/TR	Language/Art. No. RUS/EST/LT/LV/ GR/H/SLO	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ⁺⁺ mmol/L	Mg ⁺⁺ mmol/L	Cl ⁻ mmol/L	HCO ₃ ⁻ mmol/L	inorganic Phosphate mmol/L	Glucose g/L
Ci-Ca® Dialysate K2 <i>Plus</i>	F00001624	F00003579	133	2	0	1	115.75	20	1.25	1
Ci-Ca® Dialysate K4 <i>Plus</i>	F00001625	F00003580	133	4	0	1	117.75	20	1.25	1

- Sodium and bicarbonate are reduced to compensate for the systemic infusion of sodium citrate
- Calcium-free to minimise the citrate requirement
- Slightly increased magnesium as citrate also chelates with magnesium

Backed by Experts

- KDIGO AKI Guideline ², 2012: 'We suggest using CRRT, rather than standard intermittent RRT, for hemodynamically unstable patients.'
- Oudemans-van Straaten *et al.* ³, 2011: 'The full advantage of citrate anticoagulation can only be realized, however, if its risks are well controlled by means of a sound protocol, well-trained staff, and the incorporation of a citrate module in the CRRT device.'
- Zhang & Hongying ⁴, 2012: 'Regional citrate anticoagulation is superior to systemic heparin concerning circuit life and risk of bleeding.'
- Link *et al.* ⁶, 2012: 'Regional citrate anticoagulation is an effective and safe alternative to heparin.'
- Morgera *et al.* ⁷, 2009: '[The Ci-Ca CVVHD protocol] enabled an effective treatment of acute renal failure and excellent control of acid-base status as well on the systemic ionized calcium in combination with negligible clotting issues.'

Backed by Experts

- Kalb *et al.* ⁸, 2013: 'In conclusion, the [Ci-Ca CVVHD] protocol is effective to guarantee very long filter running times enabling the delivery of high dialysis doses.'
- Schultheiß *et al.* ⁹, 2012: 'The sodium balance was stable during [Ci-Ca] CVVHD treatment, with sodium values being within the reference range [...] in 91 % of runs after 72 hours.'
- Schmidt *et al.* ¹⁰, 2012: 'The EMiC2 provides superior middle molecule elimination compared to the AV 1000S dialyzer without significant difference in albumin elimination.'
- Raimundo *et al.* ¹¹, 2013: 'In conclusion, [the Ci-Ca CVVHD] protocol with target systemic ionised calcium concentration in the physiologic range maintains stable PTH levels.'
- Rimmele *et al.* ¹², 2012: 'The removal of middle molecular weight molecules is higher with SHF-HD [i.e. CVVHD using EMiC2]. Albumin loss was limited in both groups, even with SHF-HD.'

Benefits of Regional Citrate Anticoagulation

- **Delivery of prescribed dialysis dose** – through long and reliable filter running times ⁸
- **Continuous uninterrupted treatment** – through minimised clotting events and fewer treatment interruptions ^{4, 13, 18}
- **Reduced risk of bleeding** - compared to systemic anticoagulation with heparin, allowing treatment even where heparin is contraindicated e.g. HIT ^{4, 15, 16, 17, 18}
- **Fewer blood transfusions** - related to less bleeding incidents and re-priming of circuits) compared to heparin ^{15, 16, 17, 18}
- **Nursing time released for patient care** – through reduction in machine interactions when compared with heparin ^{8, 19}
- **Prolonged lifespan of the extracorporeal circuit reducing circuit usage** – compared to heparin ^{4, 13, 18}
- **Lower blood flow rates with CVVHD compared to CVVH** - smaller access catheter to be used

Patient Selection

All patients .. but especially

- Patient's pre- or postoperatively where systemic anticoagulation may be contraindicated
- Patients with significant coagulopathy – related to sepsis, large volume transfusion etc.
- Trauma patients with potentially undiagnosed bleeding points
- Patients where surgical wound healing may be compromised by systemic anticoagulation
- Patients with profound cardiovascular instability – for whom high blood flows would be detrimental

Regional Citrate Anticoagulation in CRRT


Management of Citrate Anticoagulation with multiFiltrate Ci-Ca[®]



Fully Integrated Citrate & Calcium Pumps to Safeguard



Integrated Citrate & Calcium Lines to Safeguard



Mechanically Different Citrate & Calcium Connectors to Safeguard



Citrate connector



Calcium connector

Colour Coded Connectors to Safeguard




Maintaining Dialysis and Blood Flow Rates

Overall dose \approx Dialysate flow	Blood flow	Citrate dose	Calcium dose
1600 mL/h	80 mL/min	4.0 mmol/L	1.7 mmol/L
2000 mL/h	100 mL/min	4.0 mmol/L	1.7 mmol/L
2200 mL/h	110 mL/min	4.0 mmol/L	1.7 mmol/L
2600 mL/h	130 mL/min	4.0 mmol/L	1.7 mmol/L
3000 mL/h	150 mL/min	4.0 mmol/L	1.7 mmol/L
3200 mL/h	160 mL/min	4.0 mmol/L	1.7 mmol/L
3600 mL/h	180 mL/min	4.0 mmol/L	1.7 mmol/L

- The dialysate to blood flow ratio should be at a constant numerical ratio of 20:1 (corresponding to a physical flow ratio of 1.3)
- If the patient's acid-base balance is deranged this ratio can be adjusted

Display During Citrate Anticoagulation



Post-Filter Ionised Calcium

Postfilter <i>ionised</i> calcium [mmol/L]	Change of the citrate dose (citrate/blood)
> 0.40	Increase by 0.2 mmol/L and inform physician
0.35–0.40	Increase by 0.1 mmol/L
0.25–0.34	No change (typical target range)
0.20–0.24	Decrease by 0.1 mmol/L
< 0.20	Decrease by 0.2 mmol/L and inform physician

Adjust the Citrate Flow

CVVHD Ci-Ca Treatment time: 0:00 h:min
Balance: 0.00 l **Off**
Cont. heparin adm.: **Off**

Pressure arterial	Pressure venous	TMP	Dialysate	Blood flow	Ultrafiltration
-20	100	20	2000 ml/h	100 ml/min	0 ml/h
			Citrate/blood		
			Citrate Flow 176ml/h 4.0mmol/l		
			Calcium/filtrate		
			Ca-flow 39ml/h 1.7 mmol/l		
				Cont heparin adm. Off	
				Heparin bolus Bolus	

Pressure arterial: mmHg scale from -200 to 300, value -20.
Pressure venous: mmHg scale from 0 to 500, value 100.
TMP: mmHg scale from 0 to 500, value 20.

20

Treatment modes | Preparation | **Treatment parameters** | Treatment | End of treatment | System parameters

Systemic Ionised Calcium

Systemic <i>ionised</i> calcium [mmol/L]	Change of the calcium dose (calcium/filtrate)
> 1.35	Decrease by 0.4 mmol/L and inform physician
1.21 – 1.35	Decrease by 0.2 mmol/L
1.12 – 1.20	No change (typical target range)
1.00 – 1.11	Increase by 0.2 mmol/L
< 1.00	Increase by 0.4 mmol/L and inform physician

Adjust the Calcium Flow

CVVHD Ci-Ca Treatment time: 0:00 h:min
Balance: 0.00 l **Off**
Cont. heparin adm.: **Off**

Pressure arterial	Pressure venous	TMP	Dialysate	Blood flow
-20 mmHg	100 mmHg	20 mmHg	2000 ml/h	100 ml/min
			Citrate/blood	Ultrafiltration
			Citrate Flow 176ml/h 4.0 mmol/l	0 ml/h
			Calcium/filtrate	Cont heparin adm.
			Ca-flow 39ml/h 1.7 mmol/l	Off
				Heparin bolus
				Bolus

mmHg 300 200 100 0 -100 -200
mmHg 500 400 300 200 100 0
mmHg 500 400 300 200 100 0

20

Treatment modes Preparation Treatment parameters Treatment End of treatment System parameters

Venous bubble catcher level

↓ ↑

Regional Citrate Anticoagulation in CRRT

Clinical Signs of Poor Citrate Metabolism and the Management Strategy With multiFiltrate Ci-Ca®



Potential Complications of Citrate Anticoagulation

- Heparin Induced Thrombocytopenia occurs in 3 – 5% of patients receiving standard heparin therapy ²⁰
- Reduced citrate metabolism occurs in 3.6% of patients ²¹


Poor Citrate Metabolism

The clinically diagnosable parameters that indicate impaired citrate metabolism are:

- Decreasing systemic ionised calcium level despite increasing calcium supplementation as per protocol
- Increased systemic total calcium
- A total calcium to ionised calcium ratio above 2.25 (this will be checked by the renal team / consultant)
- Increasing metabolic acidosis


If the patient demonstrates any of the above clinical signs, or the calcium infusion rate is running at greater than 2.1mmol/L, this could be indicative of poor citrate metabolism

Possible Strategy for Poor Citrate Metabolism



Revised Sliding Scale for Poor Citrate Metabolism

Postfilter <i>ionised</i> calcium [mmol/L]	Change of the citrate dose (citrate/blood)
>0.40	Increase by 0.2 mmol/L and inform physician
0.35–0.40	Increase by 0.1 mmol/L
0.25–0.34	No change (typical target range)
0.20–0.24	Decrease by 0.1 mmol/L
<0.20	Decrease by 0.2 mmol/L and inform physician



Postfilter <i>ionised</i> calcium [mmol/L]	Change of the citrate dose (citrate/blood)
>0.40	Increase by 0.1 mmol/L
0.35–0.40	No change (typical target range)
0.25–0.34	Decrease by 0.1 mmol/L
0.20–0.24	Decrease by 0.2 mmol/L and inform physician

Regional Citrate Anticoagulation in CRRT

Using Citrate to Manipulate Acid-Base Balance



Maintaining Dialysis and Blood Flow Rates

Overall dose \approx Dialysate flow	Blood flow	Citrate dose	Calcium dose
1600 mL/h	80 mL/min	4.0 mmol/L	1.7 mmol/L
2000 mL/h	100 mL/min	4.0 mmol/L	1.7 mmol/L
2200 mL/h	110 mL/min	4.0 mmol/L	1.7 mmol/L
2600 mL/h	130 mL/min	4.0 mmol/L	1.7 mmol/L
3000 mL/h	150 mL/min	4.0 mmol/L	1.7 mmol/L
3200 mL/h	160 mL/min	4.0 mmol/L	1.7 mmol/L
3600 mL/h	180 mL/min	4.0 mmol/L	1.7 mmol/L

- The dialysate to blood flow ratio should be at a constant numerical ratio of 20:1 (corresponding to a physical flow ratio of 1.3)
- If the patient's acid-base balance is deranged this ratio can be adjusted

Managing an Acidosis

- First review the patient, treat any underlying condition and ensure dialysate dose is appropriate
- Check that the dialysate and blood flow rates are set according to the protocol

To correct an acidosis *either*:

- Decrease the dialysate flow (a decrease of 20% will increase the serum bicarbonate level by approximately 4mmol/L)

or

- Increase the blood flow rate (an increase of 20% to the blood flow rate will increase the serum bicarbonate level by approximately 4mmol/L)

Managing an Alkalosis

- First review the patient, treat any underlying condition and ensure dialysate dose is appropriate
- Check that the dialysate and blood flow rates are set according to the protocol.

To correct an alkalosis *either*:

- Increase the dialysate flow (an increase of 20% will decrease the serum bicarbonate level by approximately 4mmol/L)

or

- Decrease the blood flow rate (a decrease of 20% to the blood flow rate will decrease the serum bicarbonate level by approximately 4mmol/L)

Thank You

